


Using M-PEA to measure Rapid Light Curves

A Novel Application



Using M-PEA to measure Rapid Light Curves - a novel application

Summary

Rapid Light Curves (RLCs) have been used in photosynthesis research since the late 1990s, predominantly using Pulse Amplitude Modulation (PAM) fluorometers. This application note describes how to measure RLCs using the Hansatech Instruments M-PEA. A separate note covers RLCs using the Handy PEA+.

Light curves in photosynthesis research

Various types of light curves have been used for many years in the study of photosynthetic performance because they can assess both the plant's present photosynthetic capacity and potential activity across a range of ambient light intensities (Ralph & Gademann 2005).

$P - E$ curves (or photosynthesis – irradiance curves, Figure 1) were traditionally made by the measurement of oxygen evolution rates at different ambient light intensities or from C isotope incorporation (Houliez et al. 2017).




Figure 1. An illustration of a photosynthetic light-response curve, i.e., the dependency of the rate of photosynthesis on irradiance. The α -slope shows the maximum efficiency of photosynthesis ($P - E$). The intersect between the α -slope and maximum rate of photosynthesis P_{max} indicates the light saturation (I_K or E_K); surplus light energy over saturation is dissipated. At supra-optimum irradiance, the photosynthetic rate declines (dotted part of the light-response curve); this phenomenon is usually called down-regulation or photoinhibition (taken from Masojádek et al. 2021)

Steady-state light curves (SSLCs) are a chlorophyll fluorescence-based measurement protocol plotting the Electron Transport Rate (ETR) as a function of Photon Flux Density (PPFD) and are comparable to traditional $P - E$ curves (Houliez et al. 2017).

SSLCs take measurements of photosynthetic rates during several light steps at increasing intensities. The duration of each step is sufficiently long to achieve steady-state photosynthetic rates during that step. SSLC protocols can typically take anything from several minutes to a few hours to complete. Under field conditions, this presents significant challenges when comparing results between different plants, as varying factors such as time of day, associated diurnal changes of the plant and dynamic weather conditions must be considered (Rascher et al. 2000).

Rapid Light Curves (RLCs) can be used to provide detailed information relating to the saturation characteristics of electron transport through photosystem II (PSII) in addition to overall photosynthetic performance (Ralph & Gademann 2005).

RLCs consist of a series of relatively short (<30s, typically 10s) light steps, with the light intensity increasing at each step. Each light step is interspersed by a saturating pulse where the ETR value is calculated. Unlike the $P - E$ light response curves or the SSLC, RLC protocols do not achieve steady-state photosynthetic rates during the light steps. In contrast to $P - E$ curves, which provide an indication of optimal state of photosynthetic performance independent of light history, RLCs indicate the current state of photosynthetic performance. Because steady-state conditions are not reached in RLCs, they reflect the light-acclimation state in the period leading up to the measurement and also longer-term light history (Ralph & Gademann 2005).

With ETR plotted against PPFD, RLC show three distinct regions (Figure 2):

1. Light-limited region

Photosynthetic rates are limited by low light levels in the light-limited region. The parameter α indicates the slope of the rise of ETR vs. PPFD and is proportional to efficiency of light capture (effective quantum yield or Φ_{PSII}) (Schreiber 2004).

2. Light-saturated region

During this phase, the capacity of the electron transport chain limits the electron transport rate. The ETR vs. PPFD curve reaches a plateau where maximum ETR occurs (denoted by the parameter ETR_{max}) (Schreiber 2004). The minimum saturating irradiance, denoted by the parameter E_k (sometimes referred to as I_k), is determined by finding the intercept of α and ETR_{max} (Sakshaug et al. 1997) and can be related to quenching. Below E_k , photochemical quenching is the dominant pathway whereas non-photochemical quenching is dominant above E_k (Henley 1993).

3. Photoinhibited/Down-regulated region

In this region, where the plant is subjected to supra-saturating light intensities, the ETR vs. PPFD curve often tends to decline, which could be associated with photoinhibition (Henley 1993). This effect would be more likely to occur with traditional $P - E$ or SSLCs, where steady-state photosynthetic rates are achieved. However, as steady-state is not achieved in RLC protocols, there isn't normally sufficient time for photodamage to occur. It has been suggested that the decline of ETR at supra-saturating light intensities could be linked to dynamic down-regulation of PSII (White and Critchley 1999).




Figure 2. With ETR plotted against PPFD, RLC show three distinct regions; (A) Light-limited region, (B) Light-saturated region and (C) Photoinhibited/Down-regulated region.

Measurement of Rapid Light Curves

RLC protocols typically begin with an initial quasi-dark period of no more than 10 seconds. This allows rapid reoxidation of the primary electron acceptor (Q_A), without substantial relaxation of the non-photochemical quenching coefficient (Schreiber 2004). This contrasts with an extended, full dark-adaptation period of several minutes or more, used, for example, in measurements of Fast Chlorophyll Fluorescence Induction. Under full dark-adaptation, the light and dark reactions become inactive, leading to complete re-oxidation of the photosystems (Rascher et al. 2000). If complete re-oxidation of PSI and PSII occurred, illumination in the RLC protocol would have an induction effect, making interpretation of the data more complex.

The quasi-dark period is followed immediately by a saturating pulse. This is followed by a stepped sequence of actinic periods of increasing intensity, each followed by a saturating pulse. A typical duration of the actinic periods is 10 seconds each.




Figure 3. The M-PEA (Multi-Function Plant Efficiency Analyser) chlorophyll fluorometer.

The ideal settings for the intensities of the actinic steps will generate several data points in each of the three regions described previously. For samples that photosynthesize more efficiently under lower light conditions, it is more beneficial to measure more data points during the light-limited region for a better estimation of the α parameter. If the sample is high-light tolerant, then more data points in the light-saturated and photoinhibited regions will provide more interesting results.

Measurement of RLCs using M-PEA

The M-PEA (Figure 3) software protocol editor, and specifically the pre-acquisition feature (Figure 4), can be used to measure parameters and protocols more typically associated with a modulated fluorometer, such as Φ PSII (or Y(II)), ETR, NPQ and quenching analysis models etc.

You can start to make RLC measurements at once using the sample protocol available from the Hansatech Instruments website (Support > Applications). Alternatively, you can install and modify the sample protocol, or create a protocol from scratch.

There is an important caveat about nomenclature. Because the M-PEA control and analysis software currently assumes that a dark-adapted protocol has been performed, it will by default apply dark-adapted nomenclature to measured parameters in its analysis tools.

Hansatech Instruments have provided a Microsoft Excel® template to perform the graphical and statistical analysis required to interpret RLC data. The nomenclature can be transposed, as per the information in Table 1., after exporting the M-PEA data into the spreadsheet.

Using the sample protocol

1. Download the sample RLC protocol to your PC
2. Connect your M-PEA to the PC and switch the M-PEA power on
3. Open the M-PEA software on your PC
4. If the M-PEA is not automatically located, click the torch icon on the toolbar to search for connected instruments
5. Open the Protocol Editor from the M-PEA > Protocol Editor menu option
6. In the Protocol Editor, select File > Open from the menu
7. Navigate to the sample protocol file on your PC and click Open
8. On the left side of the Protocol Editor, click the Light Table button
9. In the Light Table window (Figure 4), the Illumination settings in the Pre-Acquisition area define the actinic light step intensities for the RLC. As discussed previously, the ideal settings for the intensities of these actinic steps will generate several data points in each of the 3 regions described previously. For samples that photosynthesize more efficiently under lower light conditions, it is more beneficial to define more steps at lower light intensities for a better estimation of the α parameter. If the sample is high-light tolerant, then more actinic steps at higher intensities will provide more interesting results
10. On the left side of the Protocol Editor window, click the Upload button

11. The software will generate a notification once the upload is complete
12. Follow the instructions in the “Running the RLC Protocol” section below.

Pre Acquisition				Acquisition				
#	Ch 2	Duration(sec)	Illum (uMoles)	Duration	Points	Illum (uMoles)		
1	<input type="checkbox"/>	<input checked="" type="checkbox"/>	10	0	<input checked="" type="checkbox"/>	0.80s	116	3500
2	<input type="checkbox"/>	<input checked="" type="checkbox"/>	10	10	<input checked="" type="checkbox"/>	0.80s	116	3500
3	<input type="checkbox"/>	<input checked="" type="checkbox"/>	10	25	<input checked="" type="checkbox"/>	0.80s	116	3500
4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	10	50	<input checked="" type="checkbox"/>	0.80s	116	3500
5	<input type="checkbox"/>	<input checked="" type="checkbox"/>	10	100	<input checked="" type="checkbox"/>	0.80s	116	3500
6	<input type="checkbox"/>	<input checked="" type="checkbox"/>	10	175	<input checked="" type="checkbox"/>	0.80s	116	3500
7	<input type="checkbox"/>	<input checked="" type="checkbox"/>	10	250	<input checked="" type="checkbox"/>	0.80s	116	3500
8	<input type="checkbox"/>	<input checked="" type="checkbox"/>	10	350	<input checked="" type="checkbox"/>	0.80s	116	3500
9	<input type="checkbox"/>	<input checked="" type="checkbox"/>	10	500	<input checked="" type="checkbox"/>	0.80s	116	3500
10	<input type="checkbox"/>	<input checked="" type="checkbox"/>	10	750	<input checked="" type="checkbox"/>	0.80s	116	3500

Figure 4. The Light Table window in the M-PEA Protocol Editor showing settings loaded from the sample protocol. The settings in the “Pre Acquisition” area of the window define the light intensities of each of the individual actinic light steps.

Defining an RLC Protocol from scratch

Please follow the steps below to define the RLC protocol in the M-PEA software.

1. From the M-PEA software, open the protocol editor
2. Click on the “Options” button and firstly input the Title of your protocol
3. Using the “Records” spin box, set the value to 10. This will generate 10 individual light steps
4. All other settings in this screen may be left unchecked (Channel 2, Far red, DF, Repeats)
5. Click on the “Light Table” button to view the individual light step settings
6. Set the “Pre Acquisition” duration setting to a period <30 seconds (10 seconds recommended)
7. The quasi-dark period is set by leaving the illumination setting for number 1 as “0”
8. In the subsequent illumination boxes, define a set of actinic illumination intensities to cover all 3 regions mentioned previously. For example, a set of actinic steps such as 10, 25, 50, 100, 175, 250, 350, 500, 750, 1200, 1750 would be sufficient for a low-light/shade plant; for a high-light tolerant plant, the higher light steps in this range may not be supra-saturating and more higher intensities should be selected as opposed to the more numerous lower light intensity steps
9. Set the recording duration to 1 second using the spin box controls

10. Set all Acquisition illumination settings to $>3500 \mu\text{mol m}^{-2} \text{s}^{-1}$
11. Upload the protocol to the M-PEA
12. Follow the instructions in the “Running the RLC Protocol” section below.

Running the RLC Protocol

1. Place an M-PEA dark-adaptation leafclip on the required sample but do not close the shutter since complete re-oxidation of the photosystems will cause the RLC to have an induction effect, making data interpretation more complex (Rascher et al. 2000).
2. Fit the M-PEA sensor to the leafclip and, as quickly as possible to avoid prolonged dark-adaptation of the sample, begin the measurement by clicking the “GO” button on the software toolbar, ensuring that the sample stays securely in position within the leafclip throughout the duration of the experiment.
3. Once complete, data will automatically be transferred to the M-PEA and the Summary tab will display information relating to the measurements made.

Analysing the RLC Data

As described previously, M-PEA software by default assumes that full dark-adaptation has been performed and applies nomenclature on that basis. The data captured during the RLC protocol is from a light-adapted sample. Therefore, for each of the measurements made during the protocol, the parameters displayed in the M-PEA software should be re-labelled as shown in Table 1. below.

Default nomenclature in M-PEA software	Nomenclature for parameters measured in RLC protocol
Fo	F
Fm	Fm'
Fv	Fq'
Fv/Fm	ΦPSII

Table 1. Transposition of default nomenclature for RLC parameters made with M-PEA

M-PEA software does not currently offer the graphical and statistical analysis required to interpret RLC data. However, Hansatech Instruments have provided a Microsoft Excel® template into which M-PEA data can be transposed.

The RLC template models a curve using the equation of Platt et al. (1980) from raw data collected with the M-PEA. The template requires the Solver Add-in to be enabled within Excel® (Please [view this link](#) for further information on how to enable Solver). The Solver add-in is used to calculate the best solution/fit using the GRG nonlinear method minimising the Sum of the Squared estimate of Errors (SSE).

The template includes some sample RLC analysis worksheets for reference (Figure 5).




Figure 5. The Rapid Light Curve Excel® template allows data from the M-PEA to be analysed by plotting Electron Transport Rate over PPFD. The Solver add-in is used to calculate the best solution/fit using the GRG nonlinear method minimising the Sum of the Squared estimate of Errors (SSE) and calculates best fit values for ETR_s , α , β , ETR_{max} and E_k .

The template can be downloaded from the Support > Applications page of the Hansatech Instruments website. The following guidelines provide a walkthrough of transposing data from M-PEA software into the RLC Excel® template. For specific commands and functionality within M-PEA software, please refer to the M-PEA system manual.

1. Ensure all measured records from the RLC protocol are highlighted in the M-PEA Summary tab
2. Open the Parameters tab and ensure that the parameters F_o and F_m are displayed. If they are not, they can be enabled from the Tools > 'Parameter View' Parameters menu option
3. Copy the blank RLC template sheet into a new worksheet for each data set you wish to model
4. Under the "Measured parameters" heading, type the F_o value from the first recording in M-PEA software into cell C15 under the heading "F"
5. Type the F_m value from the first recording in M-PEA into cell D15 under the heading "Fm'"
6. From the Protocol tab in M-PEA take the PAR value from the Pre Acquisition "Illum(uM)" column (should be 0 on the first step) and type it into cell B15 under the "PAR" heading

7. Values for Φ_{PSII} and ETR are calculated automatically. ETR is calculated using the constants defined in cells B8 and B9. PFD_a and $Fraction_{PSII}$ are set to default values of 0.84 and 0.5 respectively (Baker, 2008) but can be over-ridden if necessary
8. Enter F, F_m' and PAR values for the remaining measurements rows 16 onwards until all Φ_{PSII} and ETR values have been calculated
9. Click the "Plot" button to run Solver. This will minimise SSE and a curve will be generated. Best fit values for ETRs, α , β , ETR_{max} and E_k will be calculated.

Hansatech Instruments

Hansatech Instruments is a British, scientific instrument company located in the UK. For over 45 years, our efforts have been concentrated towards the design & manufacture of high quality instrumentation for teaching & research in the fields of cellular respiration & photosynthesis. Our instruments are now in use in a wide range of programs in more than 100 countries throughout the world & have gained an enviable reputation for quality, reliability & excellent price/performance.

Products

Hansatech Instruments product range covers a wide range of applications in the fields of photosynthesis & cellular respiration.

We manufacture oxygen measurement systems based on Clark type polarographic oxygen sensors, chlorophyll fluorescence measurement systems for both continuous excitation & pulse-modulated measurement techniques & optical instrumentation for the measurement of sample chlorophyll content.

Support

Purchasers of Hansatech Instruments products can be assured of ongoing support & prompt & efficient attention to enquiries at all times. Customers are encouraged to register their instruments on our website which allows access to our Support Ticketing System in addition to instruments manuals & software upgrades



References

- BAKER NR. (2008). Chlorophyll fluorescence: a probe of photosynthesis in vivo. *Annual Review of Plant Biology* 59, 89–113.
- HENLEY WJ (1993). Measurement and interpretation of photosynthetic light-response curves in algae in the context of photoinhibition and diel changes. *J Phycol* 29:729–738
- HOULIEZ, E., LEFEBVRE, S., LIZON, F., SCHMITT, F.G., (2017). Rapid light curves (RLC) or non-sequential steady-state light curves (N-SSLC): which fluorescence-based light response curve methodology robustly characterizes phytoplankton photosynthetic activity and acclimation status? *Mar Biol* (2017) 164:175
- MASOJÍDEK, J.; RANGLOVÁ, K.; LAKATOS, G.E.; SILVA BENAVIDES, A.M.; TORZILLO, G. (2021). Variables Governing Photosynthesis and Growth in Microalgae Mass Cultures. *Processes*, 9, 820.
- PLATT, T., GALLEGOS, C.L., HARRISON, W.G. (1980). Photoinhibition of photosynthesis in natural assemblages of marine phytoplankton. *J. Mar. Res.*, 38: 687-401
- RALPH, P.J. & GADEMANN, R. (2005). Rapid light curves: a powerful tool to assess photosynthetic capacity. *Aquat. Bot.*, 82: 222–237.
- RASCHER, U., LIEBIG, M., LÜTTGE, U., (2000). Evaluation of instant light-response curves of chlorophyll fluorescence parameters obtained with a portable chlorophyll fluorometer on site in the field. *Plant Cell Environ.* 23, 1397–1405.
- SAKSHAUG, E., BRICAUD, A., DANDONNEAU, Y., FALKOWSKI, P.G., KIEFER, D.A., LEGENDRE, L., MOREL, A., PARSLOW, J., TAKAHSHI, M., (1997) Parameters of photosynthesis: definitions, theory and interpretation of results. *J Plankton Res* 19:1637–1670
- SCHREIBER U (2004). Pulse-amplitude-modulation (PAM) fluorometry and saturation pulse method: an overview. In: Papageorgiou G, Govindjee (eds) *Chlorophyll fluorescence: a signature of photosynthesis*, vol 19. Kluwer Academic Publishers, Dordrecht, pp 279–319
- WHITE, A.J., & CRITCHLEY, C. (1999). Rapid light curves: a new fluorescence method to assess the state of the photosynthetic apparatus. *Photosynth. Res.*, 59, 63–72.

Hansatech Instruments Ltd

Narborough Road, Pentney, King's Lynn, Norfolk PE32 1JL, UK

Tel: +44 (0)1760 338877 Fax: +44 (0)1760 337303

info@hansatech-instruments.com

www.hansatech-instruments.com

Hansatech
Instruments