



PP Systems (70021-1)

CIRAS-3 PC Utility Installation

Rev	Date	By	Description
1.0	6/21/2017	A. Rothwell	For CIRAS-3 PC Utility Version 2.0
1.1	9/6/2017	A. Rothwell	Updates. Added error details.

2017 PP Systems. All Rights Reserved

CIRAS-3 PC Utility Software

Important Note. If the CIRAS-3 PC Utility software is already loaded on your PC, you must uninstall it prior to loading new versions.

Uninstall Instructions

To uninstall the CIRAS-3 Utility on your PC:

1. On your Windows 10 computer, select the Start button located at the bottom left corner of the screen.
2. Find “CIRAS-3 Utility” in the list of programs and right click.
3. Select Uninstall.
4. Restart computer.

Installation Instructions

The CIRAS-3 PC Utility software is provided as a “70021-1 setup.exe” file and a MSI file named “Ciras3_PC_Utility_vxxx.msi”.

1. Insert the USB thumb drive supplied by PP Systems into your PC and navigate to the folder \CIRAS3_PC_Utility_SW which contains the “70021-1 setup.exe” and Ciras3_PC_Utility_Vxxx.msi” files.
2. Double click on the “70021-1 setup.exe” file to begin the installation process and follow the onscreen instructions until complete.

If the error, “This setup requires the .NET Framework version 2.0.50727. Please install the .NET Framework and run this setup again. The .NET Framework can be obtained from the web. Would you like to do this now?”, select No and refer to the 70021-1 CIRAS-3 WMDC Troubleshooting for Windows 10 Application Note to install WMDC (Windows Mobile Device Center). Once WMDC is installed, try installing the CIRAS-3 PC Utility Software again.

For Vista or Window 7 Operating Systems

Once installation completes, find the shortcut icon on the desktop named “CIRAS-3 Utility”. Right click and select Properties. Select Compatibility and enable “Run as administrator”.

For Windows 10 Operating System

Once installation completes, find the shortcut icon on the desktop named “CIRAS-3 Utility”. Double click the icon to open the CIRAS-3 PC Utility software. Refer to 70021-1 CIRAS-3 Remote Display Application Note for details on connecting the CIRAS-3 to a PC.